

# BD Molecular Swab Collection Kit

## Vaginal swab specimen clinician collection and transfer procedure


### Clinician collection procedure

1. Collect swab prior to pelvic, speculum, or bimanual exam.
2. Gently slide the swab no more than **2 inches (5 cm)** into the vagina. Do **not** use lubricants or other products containing substances such as carbomers.
3. Rotate the swab for 10 to 15 seconds.
4. Withdraw the swab without touching the skin outside the vagina.


### Swab-to-tube transfer procedure

Specimens collected using the BD Molecular Collection Swab must be transferred to the BD Molecular Swab Sample Buffer Tube immediately after collection.


To transfer the sample:



1. Unscrew the cap of the BD Molecular Swab Sample Buffer Tube, taking care not to contaminate the contents or the outside of the tube. Immediately after collection, insert the BD Molecular Collection Swab into the tube so that the score mark indicated by the black line is at the lip of the tube. Prior to breaking, Molecular Swab should not be placed at bottom of the tube or processing errors could occur.



2. Carefully break the shaft at the score mark and allow the swab to drop into the tube.




3. Tightly re-cap the tube.

This guide is not intended as a replacement for the detailed information contained in the IFU. Please consult BD product labels and package inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

BD Life Sciences, 7 Loveton Circle, Sparks, MD 21152-0999, U.S.  
800.638.8663

[bd.com](http://bd.com)

BD, the BD Logo, BD COR and BD MAX are trademarks of Becton, Dickinson and Company or its affiliates.  
© 2023 BD. All rights reserved. (BD-70987 3042-WW-0923 September 2023)



4. Label tube

### Specimen stability in BD Molecular Swab Sample Buffer Tube

Vaginal swab specimens can be stored for a total of 21 days at 2–30 °C in BD Molecular Swab Sample Buffer Tubes.

#### Approved for use with:

- BD Vaginal Panel for BD MAX™ and BD COR™ Systems
- BD CTGCTV2 for BD MAX™ and BD COR™ Systems

