

**PARKVIEW HEALTH LABORATORIES • Client Supply Requisition**  
**Fax Orders To: 260-458-5683**

Facility Name: \_\_\_\_\_

Address: \_\_\_\_\_

Phone Number: \_\_\_\_\_ Contact Person: \_\_\_\_\_

Date Needed: \_\_\_\_\_

Under CMS Stark Law, Parkview Health Laboratories cannot provide supplies that are not exclusively identified to be used for the collection, transportation, processing, or storage of laboratory specimens. Items such as alcohol swabs, band-aids, and tape cannot be provided because they can be used for other services besides collection of laboratory specimens. If you have any questions or concerns regarding the status of your order, please contact 260-266-1665, Monday–Friday, 6:00 a.m.–2:30 p.m.

**LAB ONLY:**

Date Filled: \_\_\_\_\_ Filled By: \_\_\_\_\_

| DESCRIPTION  | Quantity Ordered |    | DESCRIPTION                                  | Quantity Ordered |    |
|--|------------------|----|--|------------------|----|
|  | Bx               | Ea |  | Bx               | Ea |
| <b>Blood Collection Tubes</b> (each)   |                  |    | <b>Urine Collection Supplies</b>             |                  |    |
| Gold Top Serum Separator (SST) 5 ml  |                  |    | Sterile Urine Containers – Limit # given out |                  |    |
| Red Top Clot Activator <input type="checkbox"/> 10 ml <input type="checkbox"/> 4 ml  |                  |    | Non-Sterile Urine Containers                 |                  |    |
| Lavender Top (EDTA) 3 ml   |                  |    | Vacutainer Sterile CCMS Urine Collection Kit |                  |    |
| Blue Top (Sodium Citrate for PT/PTT) 2.7 ml  |                  |    | Pediatric Urine Bags                         |                  |    |
| DARK Green Top (Sodium Heparin)  |                  |    | 24 Hour Urine Containers                     |                  |    |
| DARK Green Top (Lithium Heparin) 4 ml  |                  |    | <b>Requisitions – Attach Copy</b>            |                  |    |
| Grey Top (Sodium Fluoride/Potassium Oxalate) 4 ml  |                  |    | Clinical Requisition                         |                  |    |
| Navy Top (No Additive) 7 ml  |                  |    | Cytology/Surgical Requisition                |                  |    |
| Navy Top (EDTA) 7 ml   |                  |    | Supply Requisition                           |                  |    |
| (Serotonin) Kit  |                  |    | Other (attach copy)                          |                  |    |
| Microtainer <input type="checkbox"/> Lavender <input type="checkbox"/> Green <input type="checkbox"/> Red                  |                  |    |  |                  |    |
| <b>Blood Collection Supplies</b>   |                  |    | <b>Miscellaneous Supplies</b>                |                  |    |
| Vacutainer Adapters (orange shields)   |                  |    | Glucola 50 gm                                |                  |    |
| Needles 21Gx1"   |                  |    | Amber Bags                                   |                  |    |
| Needles 22Gx1"   |                  |    | <b>Toner Model #: _____</b>                  |                  |    |
| <b>Transport Supplies</b>  |                  |    |  |                  |    |
| Specimen Transport Biohazard Bags  |                  |    |  |                  |    |
| Sputum Collection Kits   |                  |    |  |                  |    |
| Ova & Parasite Kit (3 vials)   |                  |    | <b>Patient Label Rolls</b>                   |                  |    |
| Hemocult Kits (developer not provided)   |                  |    |  |                  |    |
| Commode Hat for Urine Collection   |                  |    |  |                  |    |
| Fit Kits   |                  |    |  |                  |    |
| <b>Transport Media</b>   |                  |    |  |                  |    |
| Misc. Aerobic Culture Collection (white label)   |                  |    |  |                  |    |
| Anaerobic Specimen Collection (hot pink label)   |                  |    |  |                  |    |
| Stool Culture Transport (yellow label)   |                  |    |  |                  |    |
| Blood Culture Sets: <input type="checkbox"/> Aerobic <input type="checkbox"/> Anaerobic <input type="checkbox"/> Pediatric |                  |    |  |                  |    |
| DNA Probe Collection for GC/Chlamydia  |                  |    |  |                  |    |
| Herpes M4 Media (red top) (Limit of 5)   |                  |    |  |                  |    |
| Cary Blair Vial (green top)  |                  |    |  |                  |    |
| GeneXpert BioFire  |                  |    |  |                  |    |
| RSV Swabs  |                  |    |  |                  |    |